

The NeXus Data Format for muon Spectroscopy and Neutron or X-ray Scattering

Mark Könnecke

Paul Scherrer Institute
Switzerland

September 23, 2012

- Introduction to NeXus
 - Why NeXus?
 - NeXus concepts and rules
- Deriving application definitions and NeXus files for an example
- Break
- Interacting with NeXus file with HDF-5 tools
- Interacting with NeXus file with the NAPI

The Predicament of the Traveling Scientist

- A different data format wherever she goes

The Predicament of the Traveling Scientist

- A different data format wherever she goes
- Spends lots of time converting formats or writing readers

The Predicament of the Traveling Scientist

- A different data format wherever she goes
- Spends lots of time converting formats or writing readers
- Waits even longer to load data from inefficient data formats

The Predicament of the Traveling Scientist

- A different data format wherever she goes
- Spends lots of time converting formats or writing readers
- Waits even longer to load data from inefficient data formats
- DA requires N files in different formats, notes, local knowledge

The Predicament of the Traveling Scientist

- A different data format wherever she goes
- Spends lots of time converting formats or writing readers
- Waits even longer to load data from inefficient data formats
- DA requires N files in different formats, notes, local knowledge
- Cannot read her collaborators data

The Predicament of the Traveling Scientist

- A different data format wherever she goes
- Spends lots of time converting formats or writing readers
- Waits even longer to load data from inefficient data formats
- DA requires N files in different formats, notes, local knowledge
- Cannot read her collaborators data
- Has to keep extra information in yet another form

- Definition of a standard data format
 - Rules
 - Validation tools
- Promotion of NeXus
 - Documentation
 - NeXus API
 - Outreach to the scientific community

- Complete data for typical use
- Extendable, add additional data as you please
- Self describing
- Easy automatic plotting
- Platform independent, public domain, efficient
- Suitable for a wild variety of applications

- Devised from three independent proposals by Jonathan Tischler, APS, Przemek Klosowski, NIST and Mark Koennecke, ISIS, PSI in 94-96
- Improved during various NOBUGS conferences
- NeXus International Advisory Committee, NIAC, since 2003
- Since 2003 yearly meetings of the NIAC
- We already considered many issues!
- Except for one year, we never had money to develop NeXus


- 1 Physical file format and API for accessing files
- 2 Rules for storing data in files
- 3 Component and application definitions
- 4 NeXus Utilities

- Portable, self describing, extendable, public domain
- Hierarchical data format 5 (HDF-5)
- HDF-5:
 - grouping support
 - on the fly compression
 - reading/writing subsets
 - first dimension appendable
 - Public domain C, F77 access library
 - Used by: NASA, Boing, Deutsche Bank, the weathermen,
- XML for those who wish to edit their data
- For historical reasons we have support for the older version HDF-4


- NeXus-API hides complex HDF API
- Transparent access to all three supported physical file formats
- ANSI-C implementation
- Bindings: C++, F77, Java, python, IDL, SWIG
- January, 4, 2010: 1311217 files processed at PSI alone
- You do not have to use NAPI!
- You can write/read NeXus file with HDF-5 tools and APIs alone

- Files
- Groups identified by name and a classname beginning with NX
- Scientific data sets
- Attributes
- Links


McStas Coordinate System



NeXus Simple Coordinate System



Polar angle is always relative to the previous component




- Learning from imageCIF: be precise enough in axis descriptions to construct transformation matrices
- Allows to calculate absolute positions of components in the laboratory coordinate systems
- Can directly convert from a detector coordinate system to vectors in Lab coordinate system
- Calculate things like impact of primary beam on detector, SAS
- Allows arbitray axis to be expressed
- Intuitively describe an instrument with angles and translations and still be able to recover absolute coordinates
- Full mapping between imageCIF and NeXus now possible

$$T = \begin{pmatrix} 1 & 0 & 0 & x \\ 0 & 1 & 0 & y \\ 0 & 0 & 1 & z \\ 0 & 0 & 0 & 1 \end{pmatrix}$$

$$T = \begin{pmatrix} 1 & 0 & 0 & x \\ 0 & 1 & 0 & y \\ 0 & 0 & 1 & z \\ 0 & 0 & 0 & 1 \end{pmatrix}$$
$$R = \begin{pmatrix} r_{11} & r_{12} & r_{13} & 0 \\ r_{21} & r_{22} & r_{23} & 0 \\ r_{31} & r_{32} & r_{33} & 0 \\ 0 & 0 & 0 & 1 \end{pmatrix}$$

Combining Transformations



- Transformations can be combined by matrix multiplications
- Individual matrices can be derived by looking at the situation when everything else is 0
- Absolute positions can be obtained by multiplying the resulting matrix with its transpose
- Defines new coordinate systems at components

type rotation or translation

offset an optional translation to apply to the axis before application

direction vector around which rotated or translated


value The angle of rotation or the length of translation

dependency The order of operations to place a component

- rotation_angle, polar_angle, rotate 0 1 0
- azimuthal_angle, rotate 0 0 1
- distance, translate 0 0 1
- chi, rotate 0 0 1
- phi rotate, 0 1 0
- NeXus polar coordinate system: rotate azimuthal_angle, rotate polar_angle, translate by distance

- Axis will be annotated with attributes:
 - `vector` The direction around which to rotate or along which to translate
 - `transformation_type` rotation or translation
 - `offset` The offset to apply before the axis
 - `offset_units` The units of offset
 - `depends_on` The name of the previous axis in the dependency chain. Or `.` when this axis is the bottom of the chain.
- Additionally: a `depends_on` field per component which gives the name of the top axis in the dependency chain

Example: A Eulerian Cradle



```
sample, NXsample
  rotation_angle
    @vector=0,1,0
    @transformation_type=rotation
  chi
    @depends_on=rotation_angle
    @vector=0,0,1
    @transformation_type=rotation
  phi
    @depends_on=chi
    @vector=0,1,0
    @transformation_type=rotation
  depends_on
    phi
```

- Special group structure which can be added to any base class
- Directly specify engineering coordinates

```
geometry:NXgeometry
  translation:NXtranslation
    translation[3]
  shape:NXshape
    shape: nxbox|nxcylinder|nxsphere
    size[]
  orientation:NXorientation
    vector[3]
```

- NeXus reserves the prefix NX for group names.
- Store as much as possible
- A NeXus file has one to many NXentry groups
- There are two types of entries: raw data and processed data
- Multiple different techniques in one file go into separate NXsubentries
- If there is only one entry, the preferred name is entry, else entry1, entry2... entryn
- If an entry or an NXsubentry conforms to an application definition, the application definitions must be stated in the entries definition field.

```
entry:NXentry
  sample:NXsample

  instrument:NXinstrument
    source:NXsource
    velocity_selector:NXvelocity_selector
    detector:NXdetector
      data[xsize,ysize], signal=1 (1)
  control:NXmonitor
    data
  data:NXdata
    link to (1)
```

```
entry:NXentry
  sample:NXsample
  processing_name:NXprocess
    program
    version
    parameters:NXparameter
      raw_file
  data:NXdata
    data[nx,ny,nz], signal=1
```



```
entry:NXentry
  sample:NXsample
  instrument:NXinstrument
  ....
  sas:NXsubentry
    sample:NXsample

    instrument:NXinstrument
      source:NXsource
      velocity_selector:NXvelocity_selector
      detector:NXdetector
        data[xsize,ysize], signal=1 (1)
    control:NXmonitor
      data
    data:NXdata
      link to (1)
```

```
entry, NXentry
  measurement: NXcollection
    positions: NXcollection
      om
      two_theta
    scalars: NXcollection
      title
      wavelength
  data: NXdata
    detector1
    mca5
```

- Supports self description and allows short names in components

- Supports self description and allows short names in components
- Name, classname pair allows for multiple components of the same type

- Supports self description and allows short names in components
- Name, classname pair allows for multiple components of the same type
- NXentry allows for multiple datasets in the same file

- Supports self description and allows short names in components
- Name, classname pair allows for multiple components of the same type
- NXentry allows for multiple datasets in the same file
- NXdata supports automatic plotting

- Supports self description and allows short names in components
- Name, classname pair allows for multiple components of the same type
- NXentry allows for multiple datasets in the same file
- NXdata supports automatic plotting
- Hierarchy offers itself naturally for organising data

- Supports self description and allows short names in components
- Name, classname pair allows for multiple components of the same type
- NXentry allows for multiple datasets in the same file
- NXdata supports automatic plotting
- Hierarchy offers itself naturally for organising data
- Take care once when writing, use n times

- Store physical values in C storage order
- Use NeXus components and dictionary names
- Missing names will be quickly accepted by the NIAC
- Names: full words separated by _
- Specify units in same format as used by UDunits
- Application definitions may restrict units

- There are situations where data has to be dumped as fast as possible in order to keep up with a high data rate. Or to save disk space.
- **Data not in C storage order:** use attributes stride and offset to describe the memory layout of the data.
- **Data needs scaling:** Use a NXformula group to specify a formula in muParser notations plus the parameters and data necessary to do the scaling.
- Details on both methods will be in the NeXus manual

Associating Axis and Data, Preferred method

```
entry:NXentry
  data:NXdata
    data[nx,ny,nz], signal=1, axes=x__axis,y__axis,z__axis
    x__axis[nx]
    y__axis[ny]
    z__axis[nz]
```

Associating Axis and Data, deprecated method

- Data and axis live in the same NXgroup

```
entry:NXentry
```

```
    data:NXdata
```

```
        data[nx,ny,nz], signal=1
```

```
        x_axis[nx], axis=1
```

```
        y_axis[ny], axis=2
```

```
        z_axis[nz], axis=3
```

```
entry:NXentry
  data:NXdata
    data[nx,ny,nz], signal=1
    x_axis[nx], axis=1, primary=1
    alternate_x_axis[nx], axis=1
    y_axis[ny], axis=2
    z_axis[nz], axis=3
```

- Preserve original dimensionality of detector, if possible
- Time-of-flight becomes last dimension
- Highly irregular detectors:

```
entry:NXentry
```

```
    instrument:NXinstrument
```

```
        detector:NXdetector
```

```
            data[ndet], signal=1
```

```
            polar_angle[ndet], axis=1
```

```
            azimuthal_angle[ndet]
```

```
            distance[ndet]
```

- Come in all shapes and sizes
- Captured by rules:
 - Store all varied parameters as arrays of length NP at the appropriate place in the NeXus hierarchy
 - For multi detectors, NP, number of scan points is always the first dimension
 - In NXdata: create links to counts and varied variables

Scan Example 1: rotating sample

```
entry:NXentry
  sample:NXsample
    rotation_angle[NP], axis=1 (1)
  instrument:NXinstrument
    detector:NXdetector
      data[NP],signal=1 (2)
  control:NXmonitor
    data[NP]
  data:NXdata
    link to (1)
    link to (2)
```


Scan Example 2: complex scan in Q

```
entry:NXentry
  sample:NXsample
    rotation_angle[NP], axis=1 (1)
    phi[NP], axis=1 (2)
    chi[NP], axis=1 (3)
    h[NP], axis=1 (4), primary=1
    k[NP], axis=1 (5)
    l[NP], axis=1 (6)
  instrument:NXinstrument
    detector:NXdetector
      data[NP], signal=1 (7)
      polar_angle[NP], signal=1 (8)
  data:NXdata
    link to (1)
    link to (2)
    link to (...)
    link to (8)
```

Scan Example 3: sample rotation, area detector

```
entry:NXentry
  sample:NXsample
    rotation_angle[NP], axis=1 (1)
  instrument:NXinstrument
    detector:NXdetector
      data[NP,xsize,ysize],signal=1 (2)
  control:NXmonitor
    data[NP]
  data:NXdata
    link to (1)
    link to (2)
```

- This is rastering a sample at different wavelengths, positions etc.
- Same treatment as scans, NP replaced by NR number of raster points
- For the common case of rastering on a 2D grid one can store [nx,ny,detdim]. Be aware, though, that this causes problems if the rasterisation is aborted in mid operation.


NeXus Component and Application Definitions


- Component definitions: dictionaries of allowed field names for the various NeXus groups
- **APPLICATION DEFINITIONS**
 - **DEFINE WHAT HAS TO BE IN A NEXUS FILE FOR A CERTAIN APPLICATION**
 - **DEFINES STANDARDS**
 - **ANOTHER VIEW: CONTRACT BETWEEN FILE PRODUCERS AND USERS ABOUT WHAT HAS TO BE IN A NEXUS FILE FOR A WELL DEFINED PURPOSE**
 - **VALIDATION BY NXVALIDATE**
- Written in NeXus Definition Language, NXDL

NXaperture	NXattenuator	NXbeam_stop
NXbeam	NXbending_magnet	NXcharacterization
NXcollimator	NXcrystal	NXdata
NXdetector	NXdisk_chopper	NXentry
NXenvironment	NXevent_data	NXfermi_chopper
NXfilter	NXflipper	NXgeometry
NXguide	NXinsertion_device	NXinstrument
NXlog	NXmirror	NXmoderator
NXmonitor	NXmonochromator	NXnote
NXorientation	NXparameters	NXpolarizer
NXprocess	NXsample	NXsensor
NXshape	NXsource	NXtranslation
NXuser	NXvelocity_selector	

- 1 Construct an application definition with advice from the NIAC
 - 2 You can also inherit from and extend an existing definition
 - 3 Cure for a year; data should be produced in the new format in this time
 - 4 After curation and review: this is the standard for this application type.
- No promises, but the NIAC may do it for you
 - Description of experiment
 - Minimum set of data items necessary form common use
 - Example data

Example: WONI





- 1 **Think!** what ought to go into the file
- 2 **Map** this into the NeXus file structure
- 3 **Cast** this mapping into a NXDL file
- 4 **Standardize** your application definition together with the NIAC

- What has to go into the file?
- Minimum data necessary for common usage scenarios
- Haggle it out with your community
- Coverage ratio: $> 80\%$ of use cases

- Common usage is Rietveld analysis or profile analysis
- Data required:
 - Title
 - Sample name
 - Wavelength
 - Counts versus two_theta
 - Monitor, for normalisation

- Consider into which NeXus group an item might belong
- Look in the base class for a suitable data field
- Link the data items required for the default plot into NXdata

entry:NXentry
 title
 definition

```
entry:NXentry
  title
  definition
  sample:NXsample
    name
```

entry:NXentry

title

definition

sample:NXsample

name

instrument:NXinstrument

monochromator:NXmonochromator

wavelength

```
entry:NXentry
  title
  definition
  sample:NXsample
    name
  instrument:NXinstrument
    monochromator:NXmonochromator
      wavelength
  detector:NXdetector
    data[ndet], signal=1 (1)
    polar_angle[ndet], axis=1 (2)
```



```
entry:NXentry
  title
  definition
  sample:NXsample
    name
  instrument:NXinstrument
    monochromator:NXmonochromator
      wavelength
    detector:NXdetector
      data[ndet], signal=1 (1)
      polar_angle[ndet], axis=1 (2)
  control:NXmonitor
    data
```

```
entry:NXentry
  title
  definition
  sample:NXsample
    name
  instrument:NXinstrument
    monochromator:NXmonochromator
      wavelength
    detector:NXdetector
      data[ndet], signal=1 (1)
      polar_angle[ndet], axis=1 (2)
  control:NXmonitor
    data
  data:NXdata
    link to (1)
    link to (2)
```

From an application definition to a real NeXus file

- The structure defined by the application definition is the minimum

From an application definition to a real NeXus file

- The structure defined by the application definition is the minimum
- Practical files strive to capture much more

From an application definition to a real NeXus file

- The structure defined by the application definition is the minimum
- Practical files strive to capture much more
- Suggested procedure:
 - Look at each of your instruments components and the matching NeXus base class
 - Add whatever you feel like adding or the instrument scientists wants to have
 - Add whatever management wants to have (may be not in a NeXus group)

From an application definition to a real NeXus file

- The structure defined by the application definition is the minimum
- Practical files strive to capture much more
- Suggested procedure:
 - Look at each of your instruments components and the matching NeXus base class
 - Add whatever you feel like adding or the instrument scientists wants to have
 - Add whatever management wants to have (may be not in a NeXus group)
- Remember: Adding more fields does not break application definition compliance!

- ONOKI: ONe Of a Kind Instrument

- ONOKI: ONe Of a Kind Instrument
- Think! what you want to store

- ONOKI: ONe Of a Kind Instrument
- **Think!** what you want to store
- **Map** the list into the appropriate NeXus base classes. It helps to look at each of the components ONOKI is assembled from and to decide what you wish to store for each of them.

- ONOKI: ONe Of a Kind Instrument
- **Think!** what you want to store
- **Map** the list into the appropriate NeXus base classes. It helps to look at each of the components ONOKI is assembled from and to decide what you wish to store for each of them.
- The next one to copy ONOKI is well advised to copy what you did NeXus file wise, otherwise she will not be able to reuse your software!

Available NeXus Application Definitions

NXARCHIVE	NXMONOPD	NXREFSCAN
NXREFTOF	NXsAS	NXSCAN
NXTAS	NXTOFRAW	NXTOMO
NXTOMOPHASE	NXxeULER	NXXKAPPA
NXXNB	NXXROT	NXIQPROC
NXTOMOPROC	NXTOFSINGLE	NXDIRECTOF
NXINDIRECTOF	NXIQPROC	NXLAUETO
NXsASTOF	NXsQOM	NXTOFRAW
NXTOFSINGLE	NXXAS	NXXASPROC

Challenge 1 in science you are supposed to do new, non standard, things. These of course cannot be easily cast into a standard.

Challenge 1 in science you are supposed to do new, non standard, things. These of course cannot be easily cast into a standard.

Challenge 2 in order to establish a standard a lot of people need to agree

- Challenge 1 in science you are supposed to do new, non standard, things. These of course cannot be easily cast into a standard.
- Challenge 2 in order to establish a standard a lot of people need to agree
- Challenge 3 a standard requires scarce scientific programming resources for adoption

Chance 1 By using a discoverable data format like NeXus, XML, HDF-5, people can at least figure out what is in the data file.

- Chance 1 By using a discoverable data format like NeXus, XML, HDF-5, people can at least figure out what is in the data file.
- Chance 2 Using predefined names from a dictionary gives meaning to the data in a file.

- Chance 1 By using a discoverable data format like NeXus, XML, HDF-5, people can at least figure out what is in the data file.
- Chance 2 Using predefined names from a dictionary gives meaning to the data in a file.
- Chance 3 Using a shared API reduces learning costs and increases application stability.

- Chance 1 By using a discoverable data format like NeXus, XML, HDF-5, people can at least figure out what is in the data file.
- Chance 2 Using predefined names from a dictionary gives meaning to the data in a file.
- Chance 3 Using a shared API reduces learning costs and increases application stability.
- Chance 4 With NeXus, HDF-5 plus professional programming techniques a DA application can read any file which contains the required data.

- Chance 1 By using a discoverable data format like NeXus, XML, HDF-5, people can at least figure out what is in the data file.
- Chance 2 Using predefined names from a dictionary gives meaning to the data in a file.
- Chance 3 Using a shared API reduces learning costs and increases application stability.
- Chance 4 With NeXus, HDF-5 plus professional programming techniques a DA application can read any file which contains the required data.
- Chance 5 Storing as much data as possible increases the likelihood that the needed data is actually on file, even for unforeseen uses.

- Chance 1 By using a discoverable data format like NeXus, XML, HDF-5, people can at least figure out what is in the data file.
- Chance 2 Using predefined names from a dictionary gives meaning to the data in a file.
- Chance 3 Using a shared API reduces learning costs and increases application stability.
- Chance 4 With NeXus, HDF-5 plus professional programming techniques a DA application can read any file which contains the required data.
- Chance 5 Storing as much data as possible increases the likelihood that the needed data is actually on file, even for unforeseen uses.
- Chance 6 Application Definitions

- Soleil: 20 out of 26 instruments do NeXus, 2 mill files
- PSI-SINQ: 11 from 16 instrument on NeXus, 1.4 Mill files
- Lujan/LANL: 11 instruments, no change, 1 million files
- PSI-SLS: 2 planned,
- KEK: 10, 6 planned
- ANSTO: 7 going to 10
- ESRF: 2 beamlines, limited to NXentry, NXcollection, NXdata, moving to 4
- HZB: 3 Neutron, 1 synchrotron, 3 planned
- SNS: 14,3 in the pipeline
- DESY: 0, 11 in 2 Jahren
- Diamond: 7 NeXus only, 17 writing, moving to 18 as primary format
- ISIS: 8 using, 20 writing, planned: 20 using
- Muons: 4 instruments

- 1 Store and archive data from a wild variety of instruments
- 2 Store processed data
- 3 Store a complete workflow from raw data to publication ready data in several NXentries in one file
- 4 Store a set of related experiments in one file
- 5 Define strict and validatable standards

- New systems tend to use NeXus
- No competitor for a general purpose data format
- Planned:
 - DECTRIS will make NeXus/HDF-5 the data format for EIGER detectors
 - Collaboration with CIF
 - Better organisation for base classes
 - Enhance NXvalidate