

NeXus Code Camp

Mark Könnecke

Paul Scherrer Institute
Switzerland

October 22, 2011

- NXclonehandle
- PHDF-driver
- PyTree-API
- C++ Tree API
- Cmake
- NXdict replacement design based on NXDL
- NXvalidate
- HDF-5 1.8.*
- Fedora installer
- 64 bit dataset sizes
- All dimensions unlimited

- WWW-site from Manual
- NeXus for the impatient
- Encoding of axis dependencies
- Handling multi dimensional scans
- Sphinx discussion
- Cleanup NeXus applications
- NXdetetcor and Dectris

- Cmake (Freddy)
- Sphinx (Pete J.)
- PHDF (Mark)
- HDF-5 1.6 to 1.8 (Freddy)
- C++ Tree API (Eugen W.)
- Public NeXus Talk (Mark K.)

- NXclonehandle
 - NAPI was made threadsafe on the last code camp
 - Each thread needs a file handle: NXclonehandle
- PHDF: driver for parallel HDF-5 to write and read really fast
- PyTree API:
 - Written by Ray Osborn
 - Still needs unit tests
 - Critical for 4.3 NAPI release

- Eugen Wintersberger presents us his C++ tree API
- NXdict replacement:
 - NXdict is an additional C-API
 - Describes path of item in NeXus file via an external file
 - NXdict calls create structure as needed
 - CHANGE: use NXDL 4 structure description
 - CHANGE: review API
 - Inspiration: CDM-API
 - Seek collaboration with CDM? Affects also C++ API

- Switch from autoconf to Cmake
- NXvalidate: needs completion
- HDF-5 1.6.* to 1.8.*
- Fedora (rpm) installer
- 64 bit sizes: Freddie has to many events
- NX_UNLIMITED all over, HDF-5: yes, HDF-4, XML: No

- Short Manual (8 pages) about what NeXus is all about
- Content
 - Key Concepts
 - NeXus Benefits
 - Where do I find information?
 - How to access my data?
 - NeXus examples
 - Use cases

- Decided: extend NeXus to allow full mapping from CBF to NeXus
- Information to encode:
 - type rotation or translation: DONE!
transformation_type attribute
 - direction vector around which to rotate or along which to translate: DONE! attribute
 - value The angle of rotation or the length of translation, DONE!
 - dependency The order of operations to place a component, to be discussed!

Expressing Axis Dependency in NeXus

- Implied: use existing NeXus coordinate system
- dependson attribute pointing to depending axis
- transform field in base classes which becomes a comma separated list of the path to the transformations required to position this component
- Create a special container to hold axis dependencies, NXdependency, to collect the dependencies in one place for easy access. This is what CIF does

```
sample, NXsample  
  rotation_angle  
  chi (dependson rotation_angle)  
  phi (dependson phi)
```

```
sample, NXsample  
    rotation_angle  
    chi  
    phi  
transform = rotation_angle, chi, phi
```

```
sample, NXsample
  rotation_angle
  chi
  phi
dependency, NXdependency
  sample/chi =
    sample/rotation_angle
  sample/phi =
    sample/chi
  instrument/detector/x_translation =
    instrument/detector/distance
  instrument/detector/distance =
    instrument/detector/polar_angle
```

- Conflicting use cases:
 - Easy plotting
 - Careful data analysis

```
data, NXdata
    data[nx,ny]
        @signal=1
    x[nx]
        @axis=1
    y[ny]
        @axis=2
```

```
data, NXdata
    data[nx,ny]
        @signal=1
    x[nx,ny]
    y[nx,ny]
```



```
data, NXdata
    data[nx,ny]
        @signal=1
    x[nx,ny]
        @axis=1,2
        @label=1
    y[nx,ny]
        @axis=1,2
        @label=2
```

```
data, NXdata
    data[nx,ny]
        @signal=1
        @axes=x,y
        @axesvalue=x__scan,y__scan
    x[nx]
        @axis=1
    y[ny]
        @axis=2
    x__scan[nx,ny]
    y__scan[nx,ny]
```

```
data, NXdata  
    data[nx,ny]  
        @signal=1
```


- Dectris (Eiger, Pilatus, Mythen) going HDF-5 with NeXus conventions
- Additions to NXdetector for this kind of detector
- Programming model:
 - Dectris writes HDF-5 file with NXdetector
 - Local DAQ-system adds beamline metadata

Prioritise!

$$T = \begin{pmatrix} 1 & 0 & 0 & x \\ 0 & 1 & 0 & y \\ 0 & 0 & 1 & z \\ 0 & 0 & 0 & 1 \end{pmatrix}$$

$$T = \begin{pmatrix} 1 & 0 & 0 & x \\ 0 & 1 & 0 & y \\ 0 & 0 & 1 & z \\ 0 & 0 & 0 & 1 \end{pmatrix}$$
$$R = \begin{pmatrix} r_{11} & r_{12} & r_{13} & 0 \\ r_{21} & r_{22} & r_{23} & 0 \\ r_{31} & r_{32} & r_{33} & 0 \\ 0 & 0 & 0 & 1 \end{pmatrix}$$

Combining Transformations



- Transformations can be combined by matrix multiplications
- Individual matrices can be derived by looking at the situation when everything else is 0
- Absolute positions can be obtained by multiplying the resulting matrix with its transpose
- Defines new coordinate systems at components
- CIF contains a duplication: vector, offset scheme

- Allows to calculate absolute positions of components in the laboratory coordinate systems
- Can directly convert from a detector coordinate system to vectors in Lab coordinate system
- Calculate things like impact of primary beam on detector, SAS
- Allows arbitray axis to be expressed
- Intuitively describe an instrument with angles and translations and still be able to recover absolute coordinates

- rotation_angle, polar_angle, rotate 0 1 0
- azimuthal_angle, rotate 0 0 1
- distance, translate 0 0 1
- chi, rotate 0 0 1
- phi rotate, 0 1 0
- NeXus polar coordinate system: rotate azimuthal_angle, rotate polar_angle, translate by distance

CIF Dependency Table

axis-id	type	equipment	dependson	vector	offset
gonio_phi	rotation	goniometer	.	1,0,0,	...
det_z	translation	detector	.	0,0,-1	0 0 0
det_y	translation	detector	det_z	0,1,0	0,0,0
det_x	translation	detector	det_y	1,0,0	0,0,0