

Object Oriented NeXus Classes

Mark Könnecke

Laboratory for Development and Methods


Paul Scherrer Institut

Simple Coordinate System

- The plane perpendicular to the main rotation axis of the sample table defines the scattering plane of the instrument.
- As components are commonly positioned using angles, a polar coordinate system is used.
- When distances are required then we assume that the sample is a zero. Distances towards the source are negative, distances behind the sample, towards the detector are positive.


Polar Angle, Distance and Rotation

- The angle between the extension of the direct beam between a given component and its previous component and the projection of a third component onto the scattering plane is the polar_angle. This corresponds to longitude in a geographical coordinate system. In scattering this is synonymous with two theta or gamma in normal beam geometry.
- Birds eye view on scattering plane:



Elevation


- Standing besides the instrument:



- Elevation corresponds to latitude in geography. In neutron scattering this is often the angle ν .


Azimuthal Angle

- Again birds eye view onto the scattering plane



McStas Coordinates

- Birds eye view again:



Miscellaneous Classes

NXreflection
+h:NX_FLOAT32:[1..n] +k:NX_FLOAT32:[1..n] +l:NX_FLOAT32:[1..n]

NXmirrormaterial
+substrate_material:String +substrate_thickness:NX_FLOAT32 +substrate_roughness:NX_FLOAT32 +coating_material:String +coating_thickness:NX_FLOAT32 +coating_roughness:NX_FLOAT32 +m_value:NX_FLOAT32


NXspectrum
+wavelength:NX_FLOAT32:[1..n] +intensity:NX_FLOAT32:[1..n]

NXlog
+starttime:NXtime +offset:NX_INT32:[1..n] +value:NX_FLOAT32:[1..n]


NXuser
+name:String +affiliation:String +address:String +e-mail:String

NXarchive
+user:NXuser +instrument_reponsible:NXuser +sample:String +proposal_id: String +experiment_stardate:NXtime +title: String


NeXus Position




NeXus Shape




NeXus Stages




Source Components




Passive Beam Line Components




Active Beam Line Components



More Beamline Components




Samples




Detectors

- Concerning data handling there are the following types of detectors:
 - single detectors
 - linear detectors
 - area detectors
 - ID detectors
- Detectors have different geometries:
 - Detectors can have regular shapes: rectangles, lines etc which are best described as such
 - Some detectors (especially @ ISIS) are highly irregular: Then we need to describe each pixel.
- For each detector pixel we need to be able to deduce:
 - The scattering angle towards the previous component
 - The elevation out of the scattering plane
 - The distance to the previous component


Non-TOF Detectors




TOF-Detectors




More Detectors



Monochromatic Instruments



TOF Instruments



Conclusion

- The current NeXus classes are messy due to lack of specialization
- More classes improve clarity
- Inheritance brings better maintainability: for instance adding NXelliptical as a NXshape does not require changes downstream
- Caveats:
 - backwards compatability difficult
 - description in XML problematic