

NeXus Rules and Details


Mark Könnecke

Paul Scherrer Institute
Switzerland


May 9, 2010

- NeXus coordinate systems
- File structural rules
- Rules for storing data items
- Rules for special applications
- Special NeXus groups


McStas Coordinate System



NeXus Simple Coordinate System



Polar angle is always relative to the previous component



CIF like Coordinates (Proposed)

- Additional fields in base classes: transform, x_translation, y_translation, z_translation, aequatorial_angle
- $P_{current} = op1 * op2 * op3 * \dots * P0$
- transform is a komma separated list of axis to apply to the component to get into position.
- Additional data attributes:
 - type: translation or rotation
 - vector: rotation or translation axis
 - offset: offset to center of component
- Use this for documentation or arbitrary axis specifications

NeXus Default Coordinates (CIF like)

- polar_angle: type= rotation, vector= 0,1,0, offset= 0, 0, -distance
- azimuthal_angle: type=rotation, vector 0,0,1
- transform: azimuthal_angle, polar_angle

- Special group structure which can be added to any base class

```

geometry:NXgeometry
    translation:NXtranslation
        translation[3]
    shape:NXshape
        shape: nxbox|size[]|cylinder|nxsphere,
orientation:NXorientation
    vector[3]

```


- NeXus reserves the prefix NX for group names.
- Store as much as possible
- A NeXus file has one to many NXentry groups
- There are two types of entries: raw data and processed data
- Multiple different techniques in one file go into separate entries
- If there is only one entry, the preferred name is entry, else entry1, entry2... entryn
- If an entry conforms to an application definition, the application definitions must be stated in the entries definition field.

```
entry:NXentry
  sample:NXsample

  instrument:NXinstrument
    source:NXsource
    velocity_selector:NXvelocity_selector
    detector:NXdetector
      data[xsize,ysize], signal=1 (1)
  control:NXmonitor
    data
  data:NXdata
    link to (1)
```

```
entry:NXentry
  sample:NXsample
  processing_name:NXprocess
    program
    version
    parameters:NXparameter
      raw_file
  data:NXdata
    data[nx,ny,nz], signal=1
```

- Store physical values
- Use NeXus components and dictionary names
- Missing names will be quickly accepted by the NIAC
- Names: full words separated by _
- Specify units in same format as used by UDunits
- Application definitions may restrict units

- NeXus *STRONGLY* prefers plain data in C storage order
- Additional attributes: linearity, offset, scaling, direction, precedence
- Same meaning as in imageCIF
- linearity:
 - offset: $V_{true} = V_{raw} + offset$
 - scaling: $V_{true} = V_{raw} * scaling$
 - scaling_offset: $V_{true} = (V_{raw}/scaling) + offset$
 - sqrt_scaled: $V_{true} = (V_{raw}/scaling)^2$
 - logarithmic_scaled: $V_{true} = 10^{(V_{raw}/scaling)}$
- direction allows to select between increasing and decreasing indices
- precedence determines storage order

- Data and axis live in the same NXgroup

```
entry:NXentry
```

```
  data:NXdata
```

```
    data[nx,ny,nz], signal=1
```

```
    x_axis[nx], axis=1
```

```
    y_axis[ny], axis=2
```

```
    z_axis[nz], axis=3
```

```
entry:NXentry
  data:NXdata
    data[nx,ny,nz], signal=1
    x_axis[nx], axis=1, primary=1
    alternate_x_axis[nx], axis=1
    y_axis[ny], axis=2
    z_axis[nz], axis=3
```

```
entry:NXentry
  data:NXdata
    data[nx,ny,nz], signal=1, axes=x__axis,y__axis,z__axis
    x__axis[nx]
    y__axis[ny]
    z__axis[nz]
```


- Preserve original dimensionality of detector, if possible
- Time-of-flight becomes last dimension
- Highly irregular detectors:

```
entry:NXentry
```

```
    instrument:NXinstrument
```

```
        detector:NXdetector
```

```
            data[ndet], signal=1
```

```
            polar_angle[ndet], axis=1
```

```
            azimuthal_angle[ndet]
```

```
            distance[ndet]
```

- Come in all shapes and sizes
- Captured by rules:
 - Store all varied parameters as arrays of length NP at the appropriate place in the NeXus hierarchy
 - For multi detectors, NP, number of scan points is always the first dimension
 - In NXdata: create links to counts and varied variables

Scan Example 1: rotating sample

```
entry:NXentry
  sample:NXsample
    rotation_angle[NP], axis=1 (1)
  instrument:NXinstrument
    detector:NXdetector
      data[NP],signal=1 (2)
  control:NXmonitor
    data[NP]
  data:NXdata
    link to (1)
    link to (2)
```

Scan Example 2: complex scan in Q

```
entry:NXentry
  sample:NXsample
    rotation_angle[NP], axis=1 (1)
    phi[NP], axis=1 (2)
    chi[NP], axis=1 (3)
    h[NP], axis=1 (4), primary=1
    k[NP], axis=1 (5)
    l[NP], axis=1 (6)
  instrument:NXinstrument
    detector:NXdetector
      data[NP], signal=1 (7)
      polar_angle[NP], signal=1 (8)
  data:NXdata
    link to (1)
    link to (2)
    link to (...)
    link to (8)
```

Scan Example 3: sample rotation, area detector

```
entry:NXentry
  sample:NXsample
    rotation_angle[NP], axis=1 (1)
  instrument:NXinstrument
    detector:NXdetector
      data[NP,xsize,ysize],signal=1 (2)
  control:NXmonitor
    data[NP]
  data:NXdata
    link to (1)
    link to (2)
```

- This is rastering a sample at different wavelengths, positions etc.
- Same treatment as scans, NP replaced by NR number of raster points
- For the common case of rastering on a 2D grid one can store [nx,ny,detdim]. Be aware, though, that this causes problems if the rasterisation is aborted in mid operation.

```
name_of_logged_value:NXlog  
  time[], start  
  value[]  
  description  
  minimum_value  
  maximum_value  
  average_value
```

```
note:NXlog
  author
  type
  date
  description
  file_name
  data, NX_BINARY
```


NXaperture	NXattenuator	NXbeam_stop
NXbeam	NXbending_magnet	NXcharacterization
NXcollimator	NXcrystal	NXdata
NXdetector	NXdisk_chopper	NXentry
NXenvironment	NXevent_data	NXfermi_chopper
NXfilter	NXflipper	NXgeometry
NXguide	NXinsertion_device	NXinstrument
NXlog	NXmirror	NXmoderator
NXmonitor	NXmonochromator	NXnote
NXorientation	NXparameters	NXpolarizer
NXprocess	NXsample	NXsensor
NXshape	NXsource	NXtranslation
NXuser	NXvelocity_selector	