


NeXus Application Definitions Tutorial


Mark Könnecke

Paul Scherrer Institute
Switzerland

May 7, 2010

- Develop an application definition for WONI at INIS
- WONI = WOnderful New Instrument
- INIS = INsanely Intense Source





- A copy of the current NeXus base class definitions
- A XML editor

- 1 **Think!** what ought to go into the file
- 2 **Map** this into the NeXus file structure
- 3 **Cast** this mapping into a NXDL file
- 4 **Standardize** your application definition together with the NIAC

- What has to go into the file?
- Minimum data necessary for common usage scenarios
- Haggle it out with your community
- Coverage ratio: $> 80\%$ of use cases

- Common usage is Rietveld analysis or profile analysis
- Data required:
 - Title
 - Sample name
 - Wavelength
 - Counts versus two_theta
 - Monitor, for normalisation

- Consider into which NeXus group an item might belong
- Look in the base class for a suitable data field
- Link the data items required for the default plot into NXdata

entry:NXentry
 title
 definition

```
entry:NXentry
  title
  definition
  sample:NXsample
    name
```

entry:NXentry

title

definition

sample:NXsample

name

instrument:NXinstrument

monochromator:NXmonochromator

wavelength

```
entry:NXentry
  title
  definition
  sample:NXsample
    name
  instrument:NXinstrument
    monochromator:NXmonochromator
      wavelength
  detector:NXdetector
    data[ndet], signal=1 (1)
    polar_angle[ndet], axis=1 (2)
```

```
entry:NXentry
  title
  definition
  sample:NXsample
    name
  instrument:NXinstrument
    monochromator:NXmonochromator
      wavelength
    detector:NXdetector
      data[ndet], signal=1 (1)
      polar_angle[ndet], axis=1 (2)
  control:NXmonitor
    data
```

```
entry:NXentry
  title
  definition
  sample:NXsample
    name
  instrument:NXinstrument
    monochromator:NXmonochromator
      wavelength
    detector:NXdetector
      data[ndet], signal=1 (1)
      polar_angle[ndet], axis=1 (2)
  control:NXmonitor
    data
  data:NXdata
    link to (1)
    link to (2)
```

```
<group type="NXsource" name="'source'">  
  
</group>
```



```
<field name="data" type="NX_INT" signal="1">  
  <doc>  
    Some blabla  
  </doc>  
  <dimensions size="3">  
    <dim index="1" value="np" />  
    <dim index="2" value="number of x pixels" />  
    <dim index="3" value="number of y pixels" />  
  </dimensions>  
  <attribute name="signal" type="NX_CHAR">  
    <enumeration>  
      <item value="1" />  
    </enumeration>  
  </attribute>  
</field>
```

- There is a XML schema for NXDL, the editor helps you!

- Forward you application definition to the NIAC for review
- Correct you definition according to NIAC comments, if any
- Cure and use the definition for a year, data should be written and analysed in this year
- After a final review, this is the standard for that application. Period.

From an application definition to a real NeXus file

- The structure defined by the application definition is the minimum

From an application definition to a real NeXus file

- The structure defined by the application definition is the minimum
- Practical files strive to capture much more

From an application definition to a real NeXus file

- The structure defined by the application definition is the minimum
- Practical files strive to capture much more
- Suggested procedure:
 - Look at each of your instruments components and the matching NeXus base class
 - Add whatever you feel like adding or the instrument scientists wants to have
 - Add whatever management wants to have (may be not in a NeXus group)

From an application definition to a real NeXus file

- The structure defined by the application definition is the minimum
- Practical files strive to capture much more
- Suggested procedure:
 - Look at each of your instruments components and the matching NeXus base class
 - Add whatever you feel like adding or the instrument scientists wants to have
 - Add whatever management wants to have (may be not in a NeXus group)
- Remember: Adding more fields does not break application definition compliance!

- ONOKI: ONe Of a Kind Instrument

- ONOKI: ONe Of a Kind Instrument
- **Think!** what you want to store

- ONOKI: ONe Of a Kind Instrument
- **Think!** what you want to store
- **Map** the list into the appropriate NeXus base classes. It helps to look at each of the components ONOKI is assembled from and to decide what you wish to store for each of them.

- ONOKI: ONe Of a Kind Instrument
- **Think!** what you want to store
- **Map** the list into the appropriate NeXus base classes. It helps to look at each of the components ONOKI is assembled from and to decide what you wish to store for each of them.
- The next one to copy ONOKI is well advised to copy what you did NeXus file wise, otherwise she will not be able to reuse your software!

- Not quite ready yet
- `nxvalidate nexusfile.hdf`